



Histopathology - Tissue Processing and Paraffin Embedding CRF

Sample Transportation			
CHAIN ID (Write CHAIN ID #. Affix barcode if available)	[3] [0] [0] [0] [1] [] [] []		
Specimen Kit ID # (Write Specimen Kit ID #. Affix barcode if available)	MC- 0 _ _ _		
Date samples arrived ___ / ___ / ___ DD / MM / YYYY	Time samples arrived ___ : ___ 24 hour		
Ensure that a minimum of 4 hours has passed since collection before proceeding to tissue processing. If < 4 hours have passed, allow the formalin jars to sit in the laboratory until at least 4 hours after collection before beginning tissue processing.			
If tissue collection occurred at night, the histopathology tech will prepare the cassettes ASAP the next morning. If tissue processing cannot be conducted for a few days (i.e. weekend), the samples must be transferred to distilled water ASAP once 4 hours have passed since collection. The tissues can sit in distilled water for a few days to a few weeks until tissue processing can be conducted.			
Date samples transferred to distilled water (DW) <input type="checkbox"/> N/A	___ / ___ / ___ DD / MM / YYYY	Time samples transferred to DW <input type="checkbox"/> N/A	___ : ___ 24 hour
If samples were transferred to distilled water, place labels on the DW jars that correspond to the formalin jars the samples were transferred from.			

Tissue Processing			
Date processing begins ___ / ___ / ___ DD / MM / YYYY	Time processing begins ___ : ___ 24 hour		
Core Specimens in Formalin Jars			
Liver 05-J1 <input type="checkbox"/> N/A	Number of cores in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm		
Abdominal Organs 06-J1 <input type="checkbox"/> N/A	Number of cores in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm		
Right Thorax 07-J1 <input type="checkbox"/> N/A	Number of cores in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm		
Left Thorax 08-J1 <input type="checkbox"/> N/A	Number of cores in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm		

Brain 10-J1 <input type="checkbox"/> N/A	Number of cores in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Bone Marrow 11-J1 <input type="checkbox"/> N/A	Number of cores in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
<i>Fill in the following blanks next to "Extra Specimen #" with the corresponding label IDs affixed to the formalin jars.</i>	
Extra Specimen 1: _____ <input type="checkbox"/> Skin <input type="checkbox"/> Lymph node <input type="checkbox"/> N/A	Number of samples in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Extra Specimen 2: _____ <input type="checkbox"/> Skin <input type="checkbox"/> Lymph node <input type="checkbox"/> N/A	Number of samples in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Extra Specimen 3: _____ <input type="checkbox"/> Skin <input type="checkbox"/> Lymph node <input type="checkbox"/> N/A	Number of samples in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Extra Specimen 4: _____ <input type="checkbox"/> Skin <input type="checkbox"/> Lymph node <input type="checkbox"/> N/A	Number of samples in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm

GI Specimens in Formalin Jars	
D3/D4/jejunum 04-J1 <input type="checkbox"/> N/A (If specimens placed in distilled H2O, retain single bx in 04-J1, place MITS bxs in jar labeled with sharpie.)	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm Number of biopsies used for paraffin embedding: _____ Number of biopsies retained in formalin jar: _____
D2 05-J1 <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
D1 07-J1 <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Stomach 08-J1 <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Rectum 09-J1 <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
TI/Colon 10-J1 <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
<i>Fill in the following blanks next to "Extra Specimen #" with the corresponding label IDs affixed to the formalin jars.</i>	
Extra Specimen 1: _____ <input type="checkbox"/> Abnormal <input type="checkbox"/> Esophagus <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm

Extra Specimen 2: _____ <input type="checkbox"/> Abnormal <input type="checkbox"/> Esophagus <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm
Extra Specimen 3: _____ <input type="checkbox"/> Abnormal <input type="checkbox"/> Esophagus <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of shortest core: _____ mm
Extra Specimen 4: _____ <input type="checkbox"/> Abnormal <input type="checkbox"/> Esophagus <input type="checkbox"/> N/A	Number of biopsies in jar: _____ <input type="checkbox"/> >10 cores in jar Length of longest core: _____ mm Length of shortest core: _____ mm

Personnel Performing Processing and Embedding			
Name of histopathology technician	_____		
Name(s) of all assistant(s)	_____ Please separate names with a semi-colon _____ _____		
Date processing and embedding completed	____ / ____ / ____ DD / MM / YYYY	Time processing and embedding completed	____ : ____ 24 hour

Comments on Processing and Embedding Procedure (Optional)

Document History		
Version	Change Log	Date
V.1 Histopathology - Tissue Processing and Paraffin Embedding CRF	Original Version	27-Aug-18
V.1 Histopathology - Tissue Processing and Paraffin Embedding CRF	Added instructions to preserve extra D3/D4/jejenum biopsy in formalin jar	13-Dec-18