

Overview of Mycobacterium Growth Indicator Tubes (MGIT) including Drug Susceptibility Testing (DST)

Priya Solanki & Linzy Elton

1. Principles
2. Reagents
3. Inoculation
4. Unloading MGIT tubes
5. Processing MGIT tubes
6. Reporting results

1. Principles

- Fluorescent compound present in the silicone at the bottom of MGIT
- Fluorescent compound is sensitive to the presence of oxygen
- Initially there is large amounts of dissolved oxygen
- The dissolved oxygen decreases as the *Mycobacterium spp* grow
- This results in increased fluorescence
- MGIT incubation is at 35-37°C
- Each MGIT is monitored every 60 minutes for increasing fluorescence



Initially



Dissolved Oxygen



Little Fluorescence



Dissolved Oxygen



Increased Fluorescence



- A **positive** MGIT contains 10^5 to 10^6 colony forming units per milliliter (CFU/mL)
- A **negative** MGIT has no fluorescence/visible signs of bacterial growth at 42 days

MGIT tubes contain:

- ✓ 110 μ L fluorescent indicator
- ✓ 7 mL modified Middlebrook 7H9 broth

2. Reagents

- Growth supplement (OADC) is added to tubes to enable the rapid growth of mycobacteria

Middlebrook OADC supplement – values per L Purified Water		
Oleic Acid	0.1 g	Involved in metabolism of mycobacteria
Albumin	50.0 g	Protective agent by binding free fatty acids which may be toxic to Mycobacterium species
Dextrose	20.0 g	Energy Source
Catalase	0.03 g	Destroys toxic peroxides that may be present in media
Polyoxyethylene stearate (POES)	1.1 g	Emulsifying agent. Encourages Mycobacterium growth



2. Reagents

- Contamination is reduced when PANTA antibiotic mixture is added to the media

PANTA: per vial of lyophilised PANTA – per L purified water	
Polymyxin B	6,000 units
Amphotericin B	600 µg
Nalidixic acid	2,400 µg
Trimethoprim	600 µg
Azlocilin	600 µg

3. Inoculation

1. Reconstitute PANTA powder with growth supplement (15mL)
2. Label the MGIT with specimen number & date
3. Add 800 μ L of growth supplement/PANTA to MGIT tube
4. Add 500 μ L of concentrated specimen suspension as prepared
5. Tightly recap the tube and invert gently
6. Scan in the tube (BD BACTEC MGIT machine) – tubes automatically tested for the recommended 42 day testing protocol



All sub-culturing must take place using **bio-safety level III** practices and containment facilities


Exact practices may vary between sites due to individual requirements. This will be highlighted with a star (*)

4. Unloading MGIT tubes

- After pressing the unload positive button
- Any positive MGIT will show a flashing light in the position it has been placed in the machine
- The tube will need to be scanned out
- A positive report can be printed and annotated as required

- Negative MGIT tubes are unloaded similarly, after pressing the unload negative button

4. Unloading MGIT tubes



4. Unloading MGIT tubes

BACTEC MGIT 960 Unloaded Positives Report

Instrument Number	Current Date/Time	Temperature A	Temperature B	Temperature C	Software Version	Page Number
1	17-06-19 10:43	36.6°C	36.9°C	36.7°C	V5.02A	1

Tube Position	Accession Number	Sequence Number	Growth Unit	Tube Status	TTD	Date Positive	Protocol Length	Start of Protocol
A/H15 - H37Rv: P5	---	430187231461	2447	+	9;21	13-06-19	42	03-06-19 17:47

Must be ≥ 4 days

5. Processing tubes: Positive

1. Scan out a tube
2. Positive tubes for *Mycobacterium tuberculosis* should flag on or after 4 days
3. Positive tubes should be sub-cultured onto solid-LJ/liquid media-MGIT as required*
4. An acid-fast smear should be prepared
5. Using a sterile pipette remove approximately 100 μL from the bottom of the tube for stain preparations
6. Report preliminary results only after acid-fast smear evaluation

5. Processing tubes: Negative

1. Scan out a tube
2. Negative tubes should be inspected after the 42 day period for positivity
3. If any signs of positivity, it should be subcultured and treated as a presumptive positive, provided the acid-fast smear result is positive
4. If there is no sign of positivity, the tube should be removed from the MGIT instrument and sterilised before discarding

- Results should be reported as per the requirements of the individual institute*
- This may vary between sites*

1. Overview
2. MDR / XDR TB
3. Principle
4. Preparation of drugs
5. Preparation of MGIT
6. Using inoculum from a positive MGIT
7. Growth control preparation
8. Inoculation of MGIT
9. Loading DST
10. Interpreting DST results
11. Confirming resistant isolates
12. Error messages
13. Quality Control

Detection of resistance in *Mycobacterium tuberculosis* is of great importance for:

Effective Patient Management

Infection Control

- Drug susceptibility testing is a rapid qualitative procedure for susceptibility testing of *Mycobacterium tuberculosis*
- Treatment of TB is most commonly through a multiple drug regimen that includes: Streptomycin (STR), Isoniazid (INH), Rifampicin (RIF), Ethambutol (EMB) and Pyrazinamide (PZA)
- It is imperative the drugs prescribed for any particular patient show appropriate activity against *Mycobacterium tuberculosis*

- Resistance to the first line drugs Isoniazid and Rifampicin is increasing globally
- Resistance to these drugs defines multi drug resistant tuberculosis (MDR TB)
- MDR TB is linked to significant mortality and is a serious threat to the efficacy of TB control programs

Extensively drug resistant tuberculosis (XDR-TB) is defined as:

- MDR-TB with resistance to second line fluoroquinolones (*e.g.* moxifloxacin)
- Plus resistance to one of the second line aminoglycosides (*e.g.* kanamycin)
- XDR-TB poses a greater challenge to treatment in patients

The BD BACTEC MGIT 960 DST principle is based on growth of the *Mycobacterium tuberculosis* isolate in a drug containing tube compared to a drug free tube (Growth Control-GC)



Growth Control

Drug containing tube

Drug Susceptibility Testing Principle Overview

Bar coded tube carrier



GC Drug containing tubes
One specimen

- DST is performed using an AST (antibiotic susceptibility testing) set
- The set consists of a **drug-free growth control MGIT tube (GC)** and **one MGIT tube for each drug**, as well as a bar-coded tube carrier that holds the set
- A known concentration of drug is added to a MGIT tube, along with the specimen
- Growth of the specimen added to the drug containing tube is compared with the drug-free growth control which contains the same specimen
- The bacterial inoculum added to the drug-free growth control tube is one hundred-fold less than the inoculum added to the drug containing tubes

First line drugs:

Drug	Stock concentration of drug*	Volume added to MGIT tubes for test	Final critical concentration in MGIT tubes	Supplier (Becton Dickinson - BD)
Streptomycin (STR)	83 µg/ml	100 µl	1.0 µg/ml	BD (SIRE kit)
Isonazid (INH)	8.3 µg/ml	100 µl	0.1 µg/ml	BD (SIRE kit)
Rifampicin (RIF)	83 µg/ml	100 µl	1.0 µg/ml	BD (SIRE kit)
Ethambutol (EMB)	415 µg/ml	100 µl	5.0 µg/ml	BD (SIRE kit)

* These drugs must be reconstituted using 4 mL of sterile/deionised water to achieve concentrations indicated in the table above.

Drug	Stock concentration of drug*	Volume added to MGIT tubes for test	Final critical concentration in MGIT tubes	Supplier (Becton Dickinson - BD)
Pyrazinamide	8000 µg/ml	100 µl	100 µg/ml	BD (PZA kit)

* PZA must be reconstituted using 2.5 mL of sterile/deionised water to achieve concentrations indicated in the table above.

Second line drugs:

Drug	Stock concentration of drug*	Volume added to MGIT tubes for test	Final critical concentration in MGIT tubes	Supplier (Becton Dickinson - BD)
Moxifloxacin (MOX)	20.75 µg/ml	100 µl	0.25 µg/ml	BD (Moxifloxacin HCl)
Kanamycin (KAN)	207.5 µg/ml	100 µl	2.5 µg/ml	BD (Kanamycin Sulphate)

* MOX must be reconstituted using 3 mL of sterile/deionised water and then diluted 1:8 with sterile/deionised water to achieve concentrations indicated in the table above.

* KAN must be reconstituted using 4 mL of sterile/deionised water to achieve concentrations indicated in the table above.

Pyrazinamide (PZA) Differences

- The activity of PZA requires a lower pH of 6.0
- Green tubes provided by BD are used for PZA DST testing
- The growth control for PZA is also incubated in a green tube as the lower pH of 6.0 does not inhibit growth of *Mycobacterium tuberculosis*



1. Label as many tubes as required for testing
2. Place tubes in correct sequence
3. Aseptically add 800 μL of BACTEC MGIT Supplement (SIRE and/or PZA)
4. Aseptically pipette 100 μL of drug




800 μL supplement

100 μL drug

- Once a MGIT has flagged positive at or after 4 days, it must be used within 1-5 days
- If a MGIT is unloaded on day 0, the tube should be re-incubated for a minimum of one day (day 1)


Using inoculum from positive MGIT: CL3

Day 1 and Day 2




- Add glass beads
- Vortex 2-10 min
- Leave to settle for 30 min

Day 3, 4 and 5



Growth Control Preparation: CL3

- SIRE/MOX/KAN drug susceptibility sets require the organism suspension to be diluted 1:100 for the growth control from which 500 μL is inoculated
- PZA drug susceptibility sets require the organism suspension to be diluted **1:10** for the growth control from which 500 μL is inoculated



Organism Suspension
diluted 1:100 for GC

500 μL inoculated

Inoculation of tubes: CL3

- 500 μ L of the organism suspension is inoculated into the drug tubes
- A blood agar plate is used to check for contamination (Incubated at 35-37°C and checked at 48 hr)



500 μ L organism suspension inoculated

- SIRE/MOX/KAN have a time in protocol of 4-13 days
- PZA allows for a longer timeframe of 4-21 days



Susceptible result:

- Growth will be inhibited and fluorescence will be suppressed in the drug-containing tube
- The drug-free growth control will grow and show increasing fluorescence
- Sensitive result: Growth units are less than 100

Resistant result:

- Growth and its corresponding increase in fluorescence will be evident in both the drug-containing and the drug-free growth control tube
- Resistant result: Growth units are more than 100

Interpreting DST results

BACTEC MGIT 960
Unloaded AST Set Report

Instrument Number	Current Date/Time	Temperature			Software Version	Page Number
		A	B	C		
1	21-06-19 13:28	36.6°C	36.9°C	36.6°C	V5.02A	1

Sequence No: 439550114904 TIP: 9;19 SOP: 10-06-19 17:17 Removed Date: 21-06-19

BCC

Tube Position	Growth Unit	Status	Concentration	Drug Name
B/L15	400	C		Growth Control
B/L16	0	S	1.00 ug/mL	Streptomycin
B/L17	0	S	0.10 ug/mL	Isoniazid
B/L18	0	S	1.00 ug/mL	Rifampin
B/L19	0	S	5.00 ug/mL	Ethambutol

Sequence No: 439220066197 TIP: 9;2 SOP: 10-06-19 17:18 Removed Date: 21-06-19

BCC

Tube Position	Growth Unit	Status	Concentration	Drug Name
B/M01	400	C		Growth Control
B/M02	0	C		Undefined Drug #1 <i>MOX</i>

END OF AST SETS

Confirming Resistant Results


Sequence No: 439220065988 *BCG* TIP: 7;23 SOP: 10-06-19 17:18 Removed Date: 19-06-19

Tube Position	Growth Unit	Status	Concentration	Drug Name
B/M03	400	C	100.0 ug/mL	Growth Control
B/M04	400	R		Pyrazinamide

1. Blood agar culture (BAC) prepared - Check at 48 hr
2. Breadcrumb morphology check (no turbidity in MGIT)
3. Staining – ZN or Kinyoun (if needed)

Resistant Result Accepted	BAC has no growth
	Colony morphology is typical
Resistant Result NOT Accepted	BAC shows growth
	MGIT does not show typical morphology
	Smear shows contamination

Acid Fast Bacilli



Purple Cocci
(contamination)

Kinyoun stain showing contamination

- If AST print out shows an X
- Run has failed – GC reached 400 GU outside of acceptable timeframe
- The result is invalid and no interpretation (S/R) will be shown

- System cannot detect sufficient growth in GC tube in specified protocol time:
 1. Too little inoculum
 2. Non-viable organisms
 3. Slow growing drug-resistant strain
- On repeat sub-culture, use 3-5 days undiluted culture

- Detects indications of possible contaminated or over inoculated GC tube
- Perform BAC to check tube is pure
- On MGIT repeat sub-culture, use between 1-2 days

- Perform DST testing on H37Rv (or equivalent fully sensitive isolate)
- Each batch of reagents (drug kits and tubes) must be quality tested
- If the QC fails: all results for batch must be reviewed, new reagents purchased, and testing of samples repeated

1. Overview
2. MDR / XDR TB
3. Principle
4. Preparation of drugs
5. Preparation of MGIT
6. Using inoculum from a positive MGIT
7. Growth control preparation
8. Inoculation of MGIT
9. Loading DST
10. Interpreting DST results
11. Confirming resistant isolates
12. Error messages
13. Quality Control