

A Comparative Examination of Influenza 2 day ELISA and 3 day HA Consensus Microneutralization Assays : *A(H3N2) and A(H5N1) update*

Karen Laurie, Alan Heath, John Wood, Othmar Engelhardt

Acknowledgements

Paul-Ehrlich-Institut, Germany

Ralf Wagner
Constanze Goepfert
Nina Alex
Joanna Hammann
Britta Neumann

University of Hong Kong, Hong Kong

Malik Peiris
Mahendra Perera

Health Protection Agency, United Kingdom

Katja Hoschler
Janice Baldevarona

WHO Collaborating Centre for Reference and Research on Influenza

Chinese National Influenza Center, China

Tian Bai
Zaijiang Yu
Jianfang Zhou

Centers for Disease Control and Prevention, United States of America

Jackie Katz
Xiuhua Lu
Vic Veguilla
Feng Liu
Yaohui Bai

Molecular Epidemiology Research Division University of Siena, Italy

Emanuele Montomoli
Guilia Lapini
Sara Sbragi

Istituto Superiore di Sanita, Italy

Maria Rita Castrucci
Isabella Donatelli
Marzia Facchini

National Institute of Infectious Diseases, Japan

Noriko Kishida
Masato Tashiro
Takato Odagiri

National Institute for Biological Standards and Control, United Kingdom

Othmar Engelhardt
Diane Major
Alan Heath

John Wood, formerly NIBSC

National Institute for Health and Welfare, Finland

Thedi Ziegler

Naval Medical Research Center United States of America

Nicholas Martin
Tad Kochel

Armed Forces Health Surveillance Center United States of America

Jose Sanchez
Michael Cooper
James Cummings

Department of Virology, Norwegian Institute of Public Health, Norway

Olav Hungnes

Department of Microbiology, Faculty of Medicine, Siriraj Hospital, Mahidol University, Thailand

Pilaipan Puthavathana
Hatairat Lerdsamran
Phisanu Pooruk
Knnika Nateerom


Operational Infectious Diseases, Naval Health Research Center, United States of America

Anthony Hawksworth
Gary Brice
Ryan Ortiguerra

WHO Collaborating Centre for Reference and Research on Influenza, VIDRL, Australia

Karen Laurie
Louise Carolan



Serological Studies



- Serological studies can confirm past infection in the absence of positive virological testing and regardless of clinical presentation, thus detecting both symptomatic and asymptomatic infection.
- CONSISE Laboratory Working group formed to 'co-ordinate and standardize the international laboratory serological response'.



Microneutralization (MN) Assay

- Detects functional antibodies that block hemagglutinin binding to sialic acid residues on cells (antibodies specific to antigenic regions or stem region of hemagglutinin protein).
- Only live influenza viruses may be used → implications for viruses recommended for use only at high containment.
- Criteria for seropositivity has been established for infection with A(H1N1)pdm09 virus.
- Correlates of protection not established.
- Useful to detect antibodies specific for avian influenza viruses.
- More sensitive than the Hemagglutination Inhibition assay at lower end of titres.
- 2 day, 3 day and 7 day MN assay protocols used in different laboratories.



Microneutralization Assay Comparison

- Serological data from different locations is often compared during an outbreak to estimate the impact of a novel infection on a population.
- Standardization of assays is important for combining and comparing data.
- MN assays vary in use of protocols and determination of endpoint titres amongst laboratories worldwide.
- Comparison of MN assays using shared sera and A(H3N2) viruses found more consistency in laboratories using shorter assays, with viral antigen detection (Stephenson et al 2007).
- Knowledge on reproducibility, intra- and inter-laboratory variability of different MN assays is limited.



Summary of CONSIDE MN Assay Comparison for A(H1N1)pdm09

- Ten laboratories shared their protocols of the 2 day ELISA and 3 day hemagglutination (HA) MN assays and a consensus protocol was developed for the 3 day HA MN assay.
- Twelve laboratories from eight countries participated in the laboratory evaluation of the 2 day ELISA (WHO) and 3 day HA consensus MN assays.
- There were differences in the sensitivity of the assays between laboratories and between the MN assay methods.
- The ratio of titres between the 2 day ELISA and the 3 day HA MN assays was similar for the International Standard that was included in the study and the in-house serum samples.
- Overall, in most laboratories, there was good correlation between the results obtained using the two assay protocols.

Our results indicate that the 2 day ELISA (WHO) and 3 day HA consensus protocols for MN assays may be considered interchangeable for assays of antibodies to the influenza A(H1N1)pdm09 virus.

What about other subtypes of influenza A viruses?

Microneutralization Assay Laboratory Comparison : A(H3N2) and A(H5N1)

Study Plan

- Each laboratory should attempt to assay antibody levels in a small panel of sera using both consensus MN assay protocols: 2 day ELISA (WHO) and 3 day HA

Virus strain

- A(H3N2) - a representative wild type or reassortant virus of the vaccine strain A/Perth/16/2009 or A/Victoria/361/2011 from laboratory's own stocks.
- A(H5N1) – a representative wildtype or reassortant virus from the clade the laboratory's sera panel 'matches'.

Sera

- Approximately 10 sera comprising low, medium and high titer antibody levels.
- Sera could be from seroepidemiology studies or from vaccine studies.
- If available, the inclusion of a ferret antiserum is recommended.

Laboratory materials

- Local resources and laboratory materials shall be used
- The same cell line should be used for both 2 day and 3 day MN assays


Number of assays

- At least three comparative assays using each protocol on different days is requested

Microneutralization Assay Comparison

2 day ELISA (WHO) protocol


(Rowe et al 1999 J Clin Microbiol 37(4):937)




http://www.who.int/influenza/gisrs_laboratory/2010_12_06_serological_diagnosis_of_influenza_by_microneutralization_assay.pdf

Consensus 3 day HA protocol


1. Prepare cells to form monolayer at least 24h before required.




2. Heat inactivate sera. Dilute, add virus (100 TCID₅₀)



3. Add virus:serum to confluent MDCK monolayer. Incubate 1-2h, 37 ° C



4. Remove virus:serum, replace with media containing trypsin. Incubate 3 days @ 37 ° C, 5% CO₂



5. HA agglutination

Specific parameters were required for dilutions, calculations and incubation times. Other parameters were recommended.

A(H3N2) and A(H5N1) Microneutralization Assay Comparison – laboratories that have submitted data to date

